

IMAGING EVIDENCE




Figure 1: Original MRI scan




Figure 2: Grad-CAM activation overlay

GRAD-CAM ANALYSIS PANEL




Figure 3: (Left→Right) Original MRI | Activation Map (Jet) | GradCAM Overlay | Tumour Boundary (threshold > 0.5)
■ Red areas = high tumour probability ■ Yellow = moderate ■ Blue = low activation

TUMOUR LOCALISATION (GRAD-CAM ANALYSIS)

Parameter	Value	Notes
Estimated Area	19.3% of visible scan	Moderate
Hemisphere / Location	CENTRAL	Based on centroid x-position
Centroid Position	(289, 251) pixels	Hottest activation pixel centroid
Bounding Region	[177,111] → [409,391]	Pixel coordinates of activation box
Activation Score	70 / 100	Mean activation in region (0=low, 100=high)
GradCAM Computed	5.18s	CPU inference time

CLINICAL INTERPRETATION

The AI model identifies imaging features consistent with a glial neoplasm. The Grad-CAM localisation map highlights a region of interest with high activation intensity. Gliomas represent the most common primary brain tumours and vary widely in grade (WHO I–IV). Immediate referral to neurosurgery and neuro-oncology is strongly recommended. MRI with contrast enhancement, spectroscopy, and perfusion imaging should be obtained. Biopsy and molecular profiling (IDH status, MGMT methylation) may be required for definitive grading.

RECOMMENDATIONS

- Urgent MRI brain with gadolinium contrast enhancement
- Neurosurgical consultation within 48 hours
- Neuro-oncology multidisciplinary team review
- Biopsy / resection planning (IDH, MGMT, EGFR profiling)
- Baseline neuropsychological assessment

TECHNICAL DETAILS

Parameter	Value
Architecture	EfficientNet-B3 + ResNet50+CBAM + DenseNet121
Inference	TTA (6-view: orig, H-flip, V-flip, rot90/180/270)
Ensemble method	Confidence-weighted average + temperature scaling
GradCAM method	Standard Grad-CAM on last conv layer per backbone
Model AUC	0.9996 (macro one-vs-rest, test set n=1311)
Model Accuracy	98.78%
Sensitivity	98.68% Specificity: 99.59%
Device	CPU (no GPU detected)
Processing time	16.94s (inference) + 5.18s (GradCAM)

■ **DISCLAIMER:** This report is generated by NeuroScan AI, an AI-assisted screening tool. It is intended for informational and research purposes only and does **NOT** constitute a medical diagnosis. All findings must be reviewed and confirmed by a qualified radiologist or neurologist before any clinical decision is made. Individual model performance may vary across patient populations and imaging protocols. Reported metrics: AUC 0.9996, Accuracy 98.78%, Sensitivity 98.68%, Specificity 99.59% (test set n=1,311 images). Report generated: 2026-04-13 11:26:03 IST | Report ID: RPT-A693ADD2