

Report ID:	RPT-CC8DCC96	Generated:	2026-04-13 14:10:33 IST
Patient:	jo	Study Date:	April 13, 2026
Patient ID:	123	Modality:	MRI Brain (AI Inferred)
Model:	BrainTumorAI v2.0 Ensemble	Model AUC:	0.9996 Acc 98.78%

PRIMARY FINDING

NO
TUMOR

Diagnosis: NO TUMOR

Confidence: 77.1%

Risk Level: LOW

Processing: 16.95s (TTA 6-view ensemble)

Method: EfficientNet-B3 + ResNet50+CBAM + DenseNet121

DIFFERENTIAL DIAGNOSIS

Class	Probability	Bar Chart	Notes
Glioma	5.8%		
Meningioma	9.5%		
No Tumor	77.1%		← Primary
Pituitary	7.6%		

ENSEMBLE AGREEMENT

Model	Confidence	Bar
EfficientNet-B3	91.7%	
ResNet50+CBAM	42.9%	
DenseNet121	91.7%	
Ensemble (TTA)	77.1%	← FINAL

IMAGING EVIDENCE




Figure 1: Original MRI scan




Figure 2: Grad-CAM activation overlay

GRAD-CAM ANALYSIS PANEL




Figure 3: (Left→Right) Original MRI | Activation Map (Jet) | GradCAM Overlay | Tumour Boundary (threshold > 0.5)

■ Red areas = high tumour probability ■ Yellow = moderate ■ Blue = low activation

TUMOUR LOCALISATION (GRAD-CAM ANALYSIS)

Parameter	Value	Notes
Estimated Area	16.7% of visible scan	Moderate
Hemisphere / Location	CENTRAL	Based on centroid x-position
Centroid Position	(75, 140) pixels	Hottest activation pixel centroid
Bounding Region	[0,87] → [178,175]	Pixel coordinates of activation box
Activation Score	67 / 100	Mean activation in region (0=low, 100=high)
GradCAM Computed	4.72s	CPU inference time

CLINICAL INTERPRETATION

No significant intracranial mass lesion identified by AI analysis. The Grad-CAM map shows no focal area of high activation consistent with tumour pathology. Routine clinical follow-up is advised as clinically indicated. If symptoms persist or new neurological signs develop, repeat imaging should be considered. This result does not exclude all pathologies — clinical correlation is essential.

RECOMMENDATIONS

- 1. Correlate with clinical history and neurological examination
- 2. No urgent imaging follow-up required based on AI screening
- 3. Repeat MRI in 6–12 months if symptoms persist
- 4. Consider EEG if seizures are present

TECHNICAL DETAILS

Parameter	Value
Architecture	EfficientNet-B3 + ResNet50+CBAM + DenseNet121
Inference	TTA (6-view: orig, H-flip, V-flip, rot90/180/270)
Ensemble method	Confidence-weighted average + temperature scaling
GradCAM method	Standard Grad-CAM on last conv layer per backbone
Model AUC	0.9996 (macro one-vs-rest, test set n=1311)
Model Accuracy	98.78%
Sensitivity	98.68% Specificity: 99.59%
Device	CPU (no GPU detected)
Processing time	16.95s (inference) + 4.72s (GradCAM)

■ **DISCLAIMER:** This report is generated by NeuroScan AI, an AI-assisted screening tool. It is intended for informational and research purposes only and does **NOT** constitute a medical diagnosis. All findings must be reviewed and confirmed by a qualified radiologist or neurologist before any clinical decision is made. Individual model performance may vary across patient populations and imaging protocols. Reported metrics: AUC 0.9996, Accuracy 98.78%, Sensitivity 98.68%, Specificity 99.59% (test set n=1,311 images). Report generated: 2026-04-13 14:10:33 IST | Report ID: RPT-CC8DCC96